

**Standard FA**  
**ALL images are shot in STEREO except for**  
**begin & end of injection**

**FE-F2-Red-free**  
**SE-F2- Red-free**

**Early Phase**

**SE-F2-photo taken at “zero” timer (single image-not stereo)**

**SE-F2-photo taken at end of injection (single image-not stereo)**

**SE-F2-10-16 images @ 1-2 secs apart till about 45secs**

**Mid Phase**

**60-90 sec- SE-F2 stereo**

**SE-F1M stereo**

**2 min FE-F2 stereo**

**2-3 min SE-F2 stereo**

**Late Phase**

**5 min SE-F2 stereo**

**10 min SE-F2 stereo**

**FE-F2 stereo**

---

**FE=Fellow eye**

**SE=Study eye**

**F1M= Field 1 modified**

**F2= Field 2**

- **Film FA= shoot “reverse stereo” (right then left)**